

Blickman Adjustable Height Instrument Table | 7892SS-0



Stock Photo - Options May Vary

FACTORY NEW

Blickman 7892SS-0 - Adjustable Height Instrument Table.

Adjustable Height Instrument Tables are designed to straddle the OR table and offer twenty inches of travel to accommodate changes in table position during neuro, cardiac, and orthopedic surgery.

When used in general surgery, the table can be adjusted to accommodate the varying heights of personnel. This model can be adjusted manually from either side with a detachable crank handle.

Are you looking for a different size Instrument Table or more/less Shelves? Follow the link to view Venture Medical's catalog of Instrument Tables in a wide variety of sizes and features.

FEATURES

- Adjustable Height Instrument Table
- Designed to straddle the OR table
- Twenty inches of travel to accommodate changes in table position during neuro, cardiac, and orthopedic surgery
- When used in general surgery, the table can be adjusted to accommodate the varying heights of personnel
- This model can be adjusted manually from either side with a detachable crank handle

SPECIFICATIONS

- Model Number: 7892SS-0
- Construction: Stainless Steel
- Weight: 140 lbs
- Height: 36" to 56"
- Width: 36"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)