

Blickman Chrome 25" Round Hamper w/Casters | 2025



Stock Photo - Options May Vary

FACTORY NEW

Blickman 2025 - Chrome 25" Round Hamper w/Casters: Features include a 25" diameter ring, constructed from 1" chrome steel tubing, provides a sturdy frame, while the solid steel bag rest securely supports hamper bags. This unit rolls effortlessly on 2" casters, ensuring easy movement across floors.

Designed for long-lasting performance, it comes with a Limited Lifetime Warranty. Cover not included.

Note: Click the following link to view our Complete Catalog containing a wide variety of medical hampers, styles, sizes, chrome, stainless steel, and with or without Lids or bags. Medical Linen Hampers.

FEATURES

- All Chrome Steel Tubing
- Solid Steel Bag Rest Support
- 2" Hard Rubber Casters

SPECIFICATIONS

- Model Number: 2025
- Item Number: 096225000
- Construction: 1" Chrome Steel Tubing
- Casters: 2" Hard Rubber Casters
- Diameter: 25"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)