

Blickman Chrome Instrument/Mayo Stand | 1515



Stock Photo - Options May Vary

FACTORY NEW

Blickman 1515 - Chrome Instrument/Mayo Stand. Chrome, Single-Post Instrument/Mayo Stand features include Easily Height Adjustable with the Tru-Loc™ Friction Knob to ensure Secure Positioning of the Stainless Steel Tray. Travels easily on Hard Rubber Casters.

Blickman's Chrome Instrument/Mayo Stand comes complete with Tray and a limited Lifetime Warranty.

FEATURES

- Single-Post Chrome Mayo Stand
- Easily Adjusts
- Tru-Loc™ Friction Knob
- Stainless Steel Tray
- Travels Easily on Hard Rubber Casters

SPECIFICATIONS

- Single-Post
- Height Adjustment: 31-1/2" to 51"
- Construction: Stainless Steel Tray
- Tray Dimensions: 12-5/8"x 19-1/8"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)