

# Blickman Chrome IV Stand | 1305



Stock Photo - Options May Vary

## FACTORY NEW

Blickman 1305 - Chrome IV Stand. Value Value Value! For value-conscious customers, you will appreciate the chrome-plated construction.

This Blickman stand is offered with two hooks welded to the pole. Sits on a 4 Leg epoxy-coated base with Wall Saver™ tips to protect walls. Twist Lock Height Adjustment and stand travels on Four 2" twin-wheel hooded casters. Minimal assembly is required. Looking for IV Poles to update your existing equipment to meet your specific needs? Click the link to view our Complete Catalog of IV Poles with single and double poles, 2, 4, 6, 8 & 28 Hook/Rake poles: I.V.

Poles.

## FEATURES

- 4-Leg Base
- Epoxy-Coated Base
- Wall Saver Tips to Protect Walls
- 2 Hook Top
- Twist-Lock Height Adjustment

## SPECIFICATIONS

- 4-Leg Base
- Welded Hooks (2)
- Height Adjustment: 47 1/2" to 85 1/2"
- Hooded Twin 2" Casters
- Limited Lifetime Warranty

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)