

Blickman Decontamination Station w/SS Legs & Cross Bracing | DSHB30



Stock Photo - Options May Vary

FACTORY NEW

Blickman DSHB30 - Decontamination Station with Cross Bracing. 3 bowl soak, pre-rinse, and rinse for infection control. Blickman's new DSHB Decontamination Station w/stainless steel legs cross-bracing and stainless steel adjustable bullet feet for optimal infection control.

Blickman has the equipment you need to ensure infection control in today's high-risk healthcare environment. Blickman's Decontamination Station Features: Our durable stainless steel decontamination sinks meet the stringent AAMI (Association for the Advancement of Medical Instrumentation) guidelines for maximum sterility and utility.

These Decontamination Stations can be arranged for either left-or right-to-left workflow. Counter space is provided on both sides of the sink array.

FEATURES

- AAMI-recommended three-bowl functionality
- Ergonomic, 36" working height
- Adjustable levelers on the base of cabinets
- Unique, "no-drip" sink lip
- Integrally-welded 16-gauge cove-cornered stainless steel bowls with pitched bottoms for proper drainage

SPECIFICATIONS

- Model Number DSHB30
- Part Number 162DSHB001
- Features Unique, "no-drip" sink; base cabinets
- Dimensions 116"L x 25"D x 37"H
- Construction: Stainless Steel

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)