

Blickman Heavy Duty, 6-Leg, 4-hook, Low Center of Gravity IV Stand | 1370SS-4



Stock Photo - Options May Vary

FACTORY NEW

Blickman 1370SS-4 - Heavy Duty Low Center of Gravity IV Stand. The Blickman 1370SS-4 (6-leg, 4-hook) IV stands offer a low center of gravity base and are designed for hospital applications requiring multiple pumps or whenever a heavier base is required. The Tru-Lok friction knob secures the height position.

Comes with a Limited Lifetime Warranty. The 1370SS-4 stand with 4 hooks offers a 25-1/4" diameter low center of gravity base that offers improved stability. Heavy-duty 3" roller blade wheel clear casters for easy maneuverability.

The Blickman 1370SS-4 also features a Steel powder-coated finish base for durability. Adjusts from 52" - 90" Accommodates up to 4 infusion pumps. Color bases are available, call for details.

FEATURES

- A low center of gravity base offers improved stability
- Heavy-duty roller blade wheel clear casters for easy maneuverability
- Steel powder coated blue finish base for durability
- Color bases available (call for details)
- Adjustable Height

SPECIFICATIONS

- Base Diameter: 25-1/4"
- Construction: Steel
- Casters: Heavy-duty 3" roller blade wheel clear casters
- Finish: Powder-coated blue finish base
- Adjustable Height: 52" - 90"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)