

# Blickman Heavy Duty, 6-Leg, 6-hook, Low Center of Gravity IV Stand | 1370SS-6



Stock Photo - Options May Vary

## FACTORY NEW

Blickman 1370SS-6 - Heavy Duty Low Center of Gravity IV Stand. The Blickman 1370SS-6 (6-leg, 6-hook) IV stands offer a low center of gravity base and are designed for hospital applications requiring multiple pumps or whenever a heavier base is required. The Tru-Lok™ friction knob secures the height position.

Comes with a Limited Lifetime Warranty. The 1370SS-6 stand with 6 hooks offers a 25-1/4" diameter low center of gravity base that offers improved stability. Heavy-duty 3" roller blade wheel clear casters for easy maneuverability.

Steel powder-coated finish base for durability. The Blickman 1370SS-6 Adjusts from 52" - 90" Accommodates up to 4 infusion pumps. Note: The image may not be an exact representation.

## FEATURES

- A low center of gravity base offers improved stability
- Heavy-duty roller blade wheel clear casters for easy maneuverability
- Steel powder-coated finish base for durability
- Height Adjustable
- Accommodates up to 4 infusion pumps

## SPECIFICATIONS

- Base Diameter: 25-1/4"
- Casters: Heavy Duty 3" roller blade wheel clear casters
- Construction: Steel
- Finish: Powder Coated
- Height: 52" - 90"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)