

Blickman Heavy Duty, 6-Leg, 8-Rake, Low Center of Gravity IV Stand | 1370SS-RH



Stock Photo - Options May Vary

FACTORY NEW

Blickman 1370SS-RH - Heavy Duty Low Center of Gravity IV Stand. Blickman's new line of 1370 Six-Leg IV Stands with Low Center of Gravity Base are designed for hospital applications requiring multiple pumps or whenever a heavier base is required. The 1370SS-RH also features Eight position rake hook provides ease of use for numerous tubing lines and helps prevent entanglement.

The 1370SS-RH Features: The Tru-Loc™ friction knob secures the height position. The 25-1/4" diameter base provides stability with heavy-duty 3" casters. The Blickman 1370SS-RH has a powder-coated steel base for durability.

It adjusts from 52" to 90" and holds up to 4 pumps.

FEATURES

- 8 position rake hook provides ease of use for numerous tubing and helps prevent entanglement.
- A low center of gravity base offers improved stability
- Heavy-duty 3" roller blade wheel clear casters for easy maneuverability
- Steel powder-coated finish base for durability
- Height Adjustable

SPECIFICATIONS

- Base: Powder Coated Carbon Steel
- Casters: 3" roller blade wheel clear casters
- Finish: Steel powder coated blue finish base
- Adjustable Height: 52" - 90"
- Base Diameter: 25-1/4"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)