

Blickman Howard 7829SS Instrument Table | 30 x 16



Stock Photo - Options May Vary

FACTORY NEW

Blickman Howard 7829SS Instrument Table. Table Dimensions: 30"L x 16" W x 34"H. New and Improved Design offers complete Stainless Steel Construction. Lightweight yet extremely durable. Sanitary Closed Top-Hem prevents Bacterial Build-Up. Sealed Reinforcing Channels provide structural rigidity. 1-1/4" Diameter Legs are Fully Welded for Strength and Stability.

Are you looking for a different size Instrument Table or more/less Shelves? Follow the link to view Venture Medical's catalog of Instrument Tables in a wide variety of sizes and features.

In addition, we offer Operating Room Cabinets, Central Supply Equipment, Medical Linen Hampers, and everything in between. Click the following link to view our complete catalog of Stainless Steel Equipment.

FEATURES

- Lightweight, yet Durable
- Closed Top-Hem
- Sealed, Reinforcing Channels
- Fully Welded Legs
- 16" Clearance between Shelves

SPECIFICATIONS

- Dimensions: 30"L 16"W 34"H
- Leg Diameter: 1-1/4"
- Construction: Stainless Steel
- Limited Lifetime Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*