

Blickman Howard 7830SS Instrument Table | 20 x 16



Stock Photo - Options May Vary

FACTORY NEW

Blickman Howard 7830SS Instrument Table. Table Dimensions: 20"L x 16" W x 34" H. Design offers complete Stainless Steel Construction. Lightweight yet extremely durable. Sanitary Closed Top-Hem prevents Bacterial Build-Up.

Additional Features: The 7830SS also includes Sealed Reinforcing Channels which provide structural rigidity. Fully Welded Legs for Strength and Stability. Are you looking for a different size Instrument Table or more/less Shelves?

Follow the link to view Venture Medical's catalog of Instrument Tables in a wide variety of sizes and features.

FEATURES

- Lightweight, yet Durable
- Closed Top-Hem
- Sealed, Reinforcing Channels
- Fully Welded Legs
- 16" Clearance between Shelves

SPECIFICATIONS

- Dimensions: 20"L 16"W 34"H
- Leg Diameter: 1-1/4"
- Construction: 18-gauge (type 304) Stainless Steel
- Limited Lifetime Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*