

# Blickman Howard 7831SS Instrument Table | 33 x 18



Stock Photo - Options May Vary

## FACTORY NEW

Blickman Howard 7831SS Instrument Table. Dimensions: 33"L x 18"W x 34"H. New and Improved Design offers complete Stainless Steel Construction. Lightweight yet extremely durable.

Sanitary Closed Top-Hem prevents Bacterial Build-Up. Additional Features: The Blickman 7831SS also includes Sealed Reinforcing Channels provide structural rigidity. Legs are Fully Welded for Strength and Stability. Table travels on 2" Casters.

Are you looking for a different size Instrument Table or more/less Shelves? Follow the link to view Venture Medical's catalog of Instrument Tables in a wide variety of sizes and features.

## FEATURES

- Lightweight, yet Durable
- Closed Top-Hem
- Sealed, Reinforcing Channels
- Fully Welded Legs
- 16" Clearance between Shelves

## SPECIFICATIONS

- Dimensions: 33"L 18"w x 34"h
- Leg Diameter: 1-1/4"
- Construction: 18-gauge (type 304) Stainless Steel
- Limited Lifetime Warranty

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)