

Blickman Howard 7832SS Instrument Table | 36 x 20



Stock Photo - Options May Vary

FACTORY NEW

Blickman Howard 7832SS Instrument Table. Table Dimensions: 36"L x 20"w x 34"h. New and Improved Design offers complete Stainless Steel Construction. Lightweight yet extremely durable.

Sanitary Closed Top-Hem prevents Bacterial Build-Up. Sealed Reinforcing Channels provide structural rigidity. Legs are Fully Welded for Strength and Stability. Table travels on 2" Casters.

Are you looking for a different size Instrument Table or more/less Shelves? Follow the link to view Venture Medical's catalog of Instrument Tables in a wide variety of sizes and features.

FEATURES

- Lightweight, yet durable
- Closed Top-Hem
- Sealed, Reinforcing Channels
- Fully Welded Legs
- 16" Clearance between Shelves

SPECIFICATIONS

- Dimensions: 36"L x 20"W x 34"H
- Leg Diameter: 1-1/4"
- Construction: 18-gauge (type 304) Stainless Steel
- Limited Lifetime Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)