

Blickman Howard 7834SS Instrument Table | 48 x 20



Stock Photo - Options May Vary

FACTORY NEW

The Blickman Howard 7834SS Instrument Table features complete stainless steel construction with a lightweight yet durable design. The table measures 48"L x 20"W x 34"H and includes a sanitary closed top-hem that prevents bacterial build-up, making it ideal for medical and surgical environments.

Structural integrity is ensured through sealed reinforcing channels and fully welded legs that provide strength and stability. The table is equipped with 2" casters for easy mobility while maintaining a professional appearance suitable for clinical settings.

FEATURES

- Lightweight, yet durable
- Closed Top-Hem
- Sealed, Reinforcing Channels
- Fully Welded Legs
- 2" Casters

SPECIFICATIONS

- Dimensions: 48"L 20"W 34"H
- Leg Diameter: 1-1/4"
- Construction: 18-Gauge (type 304) Stainless Steel
- Limited Lifetime Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)