

Blickman Howard 7836SS Instrument Table | 60 x 24



Stock Photo - Options May Vary

FACTORY NEW

Blickman Howard 7836SS Instrument Table. The Blickman 7836SS Instrument Table has a new and improved design and offers 18-gauge, type 304, Stainless Steel Construction. The table dimensions are 60"l x 24"w x 34"h. Lightweight yet extremely durable. Sanitary Closed Top-Hem prevents Bacterial Build-Up.

Sealed Reinforcing Channels provide structural rigidity. 1-1/4" Diameter Legs are Fully Welded for Strength and Stability. Are you looking for a different size Instrument Table or more/less Shelves?

Follow the link to view Venture Medical's catalog of Instrument Tables in a wide variety of sizes and features.

FEATURES

- Lightweight, yet durable
- Closed Top-Hem
- Sealed, Reinforcing Channels
- Fully Welded Legs
- 16" Clearance between Shelves

SPECIFICATIONS

- Dimensions: 60"l 24"w 34"h
- Leg Diameter: 1-1/4"
- Construction: Stainless Steel
- Limited Lifetime Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)