

Blickman Howard 7837SS Instrument Table | 72 x 24



Stock Photo - Options May Vary

FACTORY NEW

Blickman Howard 7837SS Instrument Table. Table Dimensions: 72"l x 24"w x 34"h. New and Improved Design offers complete 18-gauge, type 304, Stainless Steel Construction. Lightweight yet extremely durable.

Sanitary Closed Top-Hem prevents Bacterial Build-Up. Sealed Reinforcing Channels provide structural rigidity. 1-1/4" Diameter Legs are Fully Welded for Strength and Stability. See below for more information.

Are you looking for a different size Instrument Table or more/less Shelves? Follow the link to view Venture Medical's catalog of Instrument Tables in a wide variety of sizes and features.

FEATURES

- Lightweight, yet durable
- Closed top-hem
- Sealed, reinforcing channels
- Fully welded legs
- 16" Clearance between shelves

SPECIFICATIONS

- Dimensions: 72"l x 24"w x 34"h
- Leg Diameter: 1-1/4"
- Construction: Stainless Steel
- Limited Lifetime Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)