

# Blickman Instrument / Mayo Stand | 1501



*Stock Photo - Options May Vary*

## **FACTORY NEW**

Blickman 1501 - Instrument / Mayo Stand. Features a sturdy five-leg cast aluminum base with wall savers for stability and damage protection. It includes a Tru-Loc™ friction knob for secure height adjustments and a removable stainless steel tray for easy cleaning.

The stand moves smoothly on 2-inch twin-wheel hooded casters, ensuring mobility in medical settings. Compact, durable, and backed by a limited lifetime warranty, it is a reliable tool for healthcare professionals.

Note: Click the following link to view Venture Medical's wide variety of Mayo Stand models, which best meet your needs.

## **FEATURES**

- Wall Savers to Protect Against Scuff Marks
- 5-Leg Cast Aluminum Base
- Tru-Loc™ Friction Knob Secures the Height
- Adjustable Stainless Steel tray
- Travels Easily on 2 Twin-Wheel, Hooded Casters

## **SPECIFICATIONS**

- 5-Leg Instrument Base
- Height Adjustment Range: 29 1/2" to 49"
- SS Tray Dimensions: 12-5/8" x 19-1/8"

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)