

Blickman Irrigation Stand/Tower System | IVT-4



Stock Photo - Options May Vary

FACTORY NEW

Blickman IVT-4 - Irrigation Stand/Tower System. Blickman's Revolutionary, New IV Irrigation Systems provides you with an affordable and more efficient alternative for fluid-intensive IV irrigation procedures. Additionally, the revolutionary "IVT" Irrigation Safety Systems ensure continuous fluid delivery by offering infinitely and independently adjustable heights for up to eight 3 or 5-liter IV bags.

Additional Features: They rely on gravity flow, allowing the use of standard IV tubes costing about \$5 per procedure. Their one-handed operation enables nurses to attach IV bags at comfortable heights. Additionally, bags can be raised to eight feet by pushing the Squeeze Release handles.

This ergonomic design meets OSHA standards and reduces the risk of back injuries.

FEATURES

- Affordable and more efficient alternative for fluid-intensive IV irrigation procedures
- Irrigation Safety Systems
- Independently adjustable heights for up to eight 3 or 5-liter IV bags
- Dependable gravity flow
- Approximately \$5 per procedure Easy, one-handed operation

SPECIFICATIONS

- IVT-4 - Irrigation stand/tower
- (4) channels
- Holds up to Eight 5-liter I.V. bags
- Stand: 12-3/4"W x 12-3/4"D x 77-3/4"H
- Base: 3-1/2"W x 26-3/4"D x 28-3/4"H

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)