

# Blickman Kennedy Freestanding Cabinet | 7956SS



Stock Photo - Options May Vary

## FACTORY NEW

Blickman 7956SS - Kennedy Freestanding Cabinet is a 36" wide treatment and storage cabinet with two distinct sections separated by a stainless steel work surface. The narcotic locker and three glass shelves in the upper section independently lock.

The bottom section consists of three 7" deep drawers and a closed-door compartment with one stainless steel shelf. The cabinet is mounted on four legs with stainless steel glides.

Reconfigure this cabinet with all stainless steel shelves, solid doors, a sloping top for easy cleaning, and casters.

## FEATURES

- Freestanding Cabinet
- 36" wide treatment and storage cabinet with two distinct sections separated by a stainless steel worksurface
- The narcotic locker and three glass shelves in the upper section independently lock
- The bottom section consists of three 7" deep drawers, a closed-door compartment with one stainless steel shelf
- The cabinet is mounted on four legs with stainless steel glides

## SPECIFICATIONS

- Model Number: 7956SS
- Item Number: 1527956000
- Materials: Stainless Steel, Type 300, Series, nickel-bearing, 1/4" clear tempered sheet glass.
- Steel Gauges: Sides, backs, tops, bottoms, shelf-18 gauge. Doors and drawers 20 gauge. Legs - 12 gauge.
- Construction: The cabinet consists of separate upper and lower compartments with seamless welded faces and doubled front edges.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)