

Blickman Maxi Case-Cart | CCC1-17



Stock Photo - Options May Vary

FACTORY NEW

Blickman CCC1-17 - Maxi Case-Cart. Blickman has been a leader in case-cart design for over 20 years.

Blickman introduces our new, closed-case carts with several new features to make them smooth-rolling, easy-steering, and quiet. Recommended for procedures requiring a larger volume or size of instruments and/or supplies, such as Neurosurgery and Orthopedic Surgery.

Additional Features: Choose from several of Blickman's open designs and we'll also help you design a cart to meet the unique needs of your facility, with your choice of options, such as extra handles, locks, and maintenance-free casters. All our carts are cart-washable.

FEATURES

- Maxi Case Cart has two half-size shelves for maximum flexibility
- Recommended for procedures requiring a larger volume or size of instruments and/or supplies, such as neurosurgery and orthopedic surgery
- Stainless Steel frame withstands prolonged use and has cutouts for drainage
- Body construction- 18 gauge stainless steel
- Door(s) construction - 18 gauge stainless steel

SPECIFICATIONS

- Model Number CCC1-17
- Case Cart Dimensions overall: 52"W x 29"D x 40-1/2"H
- Construction: Stainless Steel
- Body: 46"w X 25"d x 32-7/8"H
- Casters 6" diameter s/s casters, 3 with swivel; 1 with swivel lock

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)