

# Blickman Mini Case Cart | CCC3-17



Stock Photo - Options May Vary

## FACTORY NEW

The Blickman CCC3-17 Mini Case Cart is a compact, durable solution designed for labor and delivery, outpatient surgery, and instrument transfer between clinical areas. Its robust construction features an 18 gauge steel body and doors, 16 gauge top, and 11 gauge steel chassis that provide strength and stability in high-traffic hospital environments.

The cart incorporates a 14 gauge bumper channel for impact protection and added resilience. With a usable interior capacity of 9.95 cubic feet and overall dimensions of 29-5/8"W x 29"D x 40-1/2"H, this mini case cart balances compact footprint with functional storage. The design is ideal for surgical centers and hospitals where space efficiency and reliable instrument organization are essential.

## FEATURES

- Mini Case Cart
- Designed for labor and delivery, outpatient surgery
- Chassis - 11 gauge
- Bumper channel - 14 gauge
- Casters 6" diameter stainless steel: 3 Swivel, 1 with Swivel

## SPECIFICATIONS

- Model Number CCC3
- Overall Dimensions: 29-5/8"W x 29"D x 40-1/2"H
- Usable Interior: 9.95 cubic feet
- Wheel Diameter: 6"
- Estimated Shipping Weight: 175 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)