

# Blickman Motorized Straddle Adj-Height Instrument Table | 7901SS-M



Stock Photo - Options May Vary

## FACTORY NEW

Blickman 7901SS-M - Motorized Instrument Table (Liberty). The 7901SS-M is a hands-free, stainless steel motorized table for the OR. Raise or lower up to 300 lbs of instrumentation, computer monitors, and small equipment during surgery with the touch of either a foot or hand switch.

No reaching under the drape to adjust the overbed table. The 7900 Series Motorized Adjustable Instrument Tables are height adjustable from 38.75" to 58".

Are you looking for a different size Instrument Table or more/less Shelves? Follow the link to view Venture Medical's catalog of Instrument Tables in a wide variety of sizes and features.

## FEATURES

- A hands-free, stainless steel motorized table for the OR.
- Raise or lower up to 300 lbs of instrumentation, computer monitors, and small equipment during surgery
- Foot or hand switch to raise or lower
- Tables are ideal for cardiac, neuro, and orthopedic surgeries.
- 8' power cord with "Hospital Grade"

## SPECIFICATIONS

- Model: 7901SS-M
- Product Type: Motorized Operating Room Table
- Dimensions: 50.50"W x 24"D x 38.75" to 58"H
- Weight Capacity: 300 lbs
- Control: a foot or hand switch

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)