

Blickman Multi-Purpose Case Cart w/Extension Shelves | CCC2E-19



Stock Photo - Options May Vary

FACTORY NEW

Blickman CCC2E-19 - Multi-Purpose Case Cart w/Extension Shelves. The Blickman CCC2E-19 is the new and improved version of the CCC2E and CCC2E-17 models.

These Multi-Purpose Case Carts are designed with versatility in mind, offering the ideal size for most routine surgical procedures, regardless of whether your facility utilizes a container system. The cart is built for durability, with the body and doors made from 18 gauge steel, the top constructed from 16 gauge steel, and the chassis reinforced with 11 gauge steel.

Additionally, the bumper channel is crafted from 14 gauge steel, providing extra protection and stability.

FEATURES

- Multi-Purpose Case Carts
- Right size for most routine surgical procedures
- Body - 18 gauge, Door(s) - 18 gauge
- Top - 16 gauge
- Chassis - 11 gauge

SPECIFICATIONS

- Model Number CCC2E-19
- Item Number: 2293332E00
- Overall Dimensions: 42"W x 29"D x 40-1/2"H
- Construction: Body - 18 gauge, Door(s) - 18 gauge, Top - 16 gauge, Chassis - 11 gauge, and Bumper channel - 14 gauge
- Casters: 6" diameter SS

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)