

Blickman Multi-Purpose Case Cart | CCC2



Stock Photo - Options May Vary

FACTORY NEW

Blickman CCC2 – Multi-Purpose Case Cart. Multi-Purpose Case Carts are ideal for most routine surgical procedures, with or without a container system. The body uses 18-gauge steel, and the doors also use 18-gauge steel. The top features 16-gauge steel for added durability.

The chassis uses 11-gauge steel for strength and stability. The bumper channel uses 14-gauge steel for reliable protection. Heavy Duty Door Hinge improves structural integrity and allows the door to swing back 270 degrees.

Features: The Blickman CCC2E-19, an upgrade of the CCC2E and CCC2E-17, is perfect for routine surgical procedures.

FEATURES

- Multi-Purpose Case Carts
- Right size for most routine surgical procedures
- Heavy Duty Door Hinge improves structural integrity, allows the door to swing back 270 degrees
- Positive catch holds the door in an open position for cart washing
- Low maintenance, stainless steel cart

SPECIFICATIONS

- Model Number CCC2E
- Part Number 2213332596
- Overall Dimensions: 42"W x 29"D x 40 1/2"H; body: 36"W x 25"D x 32-7/8"H
- Casters 6" diameter SS
- Handles: 2 Push Handles

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)