

# Blickman Open Case Cart w/2 Wire Shelves | OCC4



Stock Photo - Options May Vary

## FACTORY NEW

Blickman OCC4 - Open Case Cart w/2 Wire Shelves. The Blickman Open Case Cart with 2 Wire Shelves offers an open architecture, making it ideal for easy cleaning and drying of stored items. This cart features two wire shelves, providing ample storage space for various supplies.

The solid bottom and bumper ensure durability, withstanding busy healthcare environments. Despite its sturdy construction, the cart remains lightweight, providing an economical solution without sacrificing reliability or strength. Its open design not only facilitates efficient cleaning but also allows for quick access to items.

This cart is perfect for use in hospitals, clinics, or surgical centers where hygiene, mobility, and storage are essential.

## FEATURES

- Open architecture
- 2 wire shelves
- Easy to clean and dry
- Lightweight yet strong
- Solid bottom

## SPECIFICATIONS

- Model Number: OCC4
- Part Number: 2410496000
- Overall Dimensions: 42"W x 29"D x 40-1/4"H
- Casters: 6" diameter stainless steel swivel, 3 swivel casters, 1 with swivel lock
- Shelves: 2

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)