

Blickman Open Case Cart w/Solid Top & Shelf | OCC3



Stock Photo - Options May Vary

FACTORY NEW

The Blickman OCC3 is an open case cart with solid top and shelf designed for efficient sterile transport in healthcare settings. Its open design facilitates easy cleaning, drying, loading, and maneuvering. The full encircling bumper and solid bottom prevent cross-contamination, making it ideal for carts that bypass non-sterile transport areas and go directly to sterile cores.

The cart features a stainless steel welded angle frame chassis that absorbs impact transmitted through wheels and bumper frames, protecting contents and preventing distortion, looseness, and rattling. A non-marking rubber bumper protects both the cart and wheels from damage during use and transport.

FEATURES

- Open Design Case Cart
- Solid Top & Shelf
- Easy to clean and dry
- Easy load and maneuver
- Full encircling bumper
- Solid bottom prevents cross-contamination
- Ideal when carts go directly to a sterile core, bypassing non-sterile transport

SPECIFICATIONS

- Model Number: OCC3
- Overall Dimensions: 42"W x 29"D x 40-1/4"H
- Casters 6" diameter, stainless steel swivel, 3 swivel casters, 1 with Swivel Lock
- Shelves: 1
- Solid Shelves
- Estimated Shipping Weight: 200 lbs
- Wheel Diameter: 5"
- Floor Space: 8.46 square feet

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)