

Blickman Space-Saver Case Cart | CCC4-19



Stock Photo - Options May Vary

FACTORY NEW

The Blickman CCC4-19 - Space-Saver Case Cart is the new and improved version of model CCC2-17. These Space-Saver Case Carts are designed to give added capacity by utilizing the cube while minimizing the floor space required by larger carts. Body - 18 gauge, Door(s) - 18 gauge, Top - 16 gauge, Chassis - 11 gauge, Bumper channel - 14 gauge.

Removable pull-out-type stainless steel wire shelves are provided on closed carts. Additional Features: Shelves slide on adjustable angle sliders, suspended from keyhole punched standards welded to the cart interior. Open carts offer either a solid stainless or wire shelf, adjustable on 4" centers.

Note: to view additional carts, click the following link [Case Carts & Operating Room Cabinets](#).

FEATURES

- Space-Saver Case Cart
- Designed to give added capacity by utilizing the cube while minimizing the floor space
- Body - 18 gauge, Door(s) - 18 gauge, Top - 16 gauge, Chassis - 11 gauge Bumper channel - 14 gauge
- Removable pull-out type stainless steel wire shelves are provided on closed carts
- Shelves slide on adjustable angle sliders, suspended from keyhole punched standards welded to the cart interior

SPECIFICATIONS

- Model Number: CCC4-19
- Part Number: 22293334000
- Dimensions: 29-5/8"W x 29"D x 55-1/2"H
- Body: 23-5/8" X 25"D X 47-7/8"H
- Casters: 6" diameter s/s, 3 swivel casters, 1 with Swivel Lock

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)