

Blickman SS Ferguson Anesthesia Utility Table w/1 Drawer | 7850SS



Stock Photo - Options May Vary

FACTORY NEW

Blickman 7850SS - Ferguson Anesthesia Utility Table w/1 Drawer. Being a popular general-use cart in healthcare environments, Blickman's superior construction includes a supporting channel underneath the top shelf to enhance its strength.

Additionally, it features a positive stop to prevent accidental drawer removal and ball-bearing slides for smooth, easy-sliding drawers. The Table features a stainless steel structure which is both durable and appropriate for a healthcare environment.

These thoughtful design elements ensure that the cart remains stable, durable, and functional, even in the most demanding healthcare settings. The cart's robust features make it a reliable choice for efficient supply organization and transport.

FEATURES

- Popular in the OR and throughout healthcare facilities
- Used as a general utility cart
- Superior construction
- Stainless Steel
- One Drawer

SPECIFICATIONS

- Model Number: 7850SS
- Part Number: 0217850000
- Dimensions: 20" L x 16" W x 34" H
- Casters Diameter: 3"
- Estimated Shipping Weight: 50lbs

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)