

Pedigo ivNow Fluid Warmer – New



Stock Photo – Options May Vary

FACTORY NEW

Pedigo ivNow Fluid Warmer – New. The ivNow1 Fluid Warmer, model P-2301-A, delivers fast and reliable warming for intravenous, injection, and irrigation solutions. It efficiently warms fluids to precise temperatures, helping clinicians provide safer, more comfortable patient care.

Each contoured warming module securely cradles 0.5-, 1-, 2-, and 3-liter solution bags, ensuring even heat distribution and consistent results. Clinicians can customize ivNow configurations by stacking or removing warming pods.

The system easily installs on countertops or mounts to walls using the included mounting plate. Each cavity operates independently with electronic temperature controls and a bright LED display, simplifying operation and monitoring.

FEATURES

- Quickly warms and maintains temperature of IV and irrigation fluids
- Accommodates 0.5-, 1-, 2-, and 3-liter solution bags
- Stackable and field-customizable pod configuration
- Usable on countertop or wall-mounted with provided plate
- Individually controlled warming cavities with LED displays

SPECIFICATIONS

- Model: P-2301-A
- Automatically warms fluid to 104°F (40°C) in 30 minutes or less
- Provides 25–57 liters of warmed fluid over a 12-hour period
- Compatible with 0.5–3 liter IV bags (per pod capacity)
- Activates automatically when IV bag is placed in pod

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)