

Pedigo ivNow2 Fluid Warmer – New



Stock Photo – Options May Vary

FACTORY NEW

Pedigo ivNow2 Fluid Warmer – New. The ivNow2 Fluid Warmer, model P-2302-A, quickly warms and maintains the temperature of intravenous, injection, and irrigation solutions before patient use. Its dual-chamber design allows two bags to warm simultaneously for greater efficiency in busy medical environments.

Each specially contoured warming cavity cradles 0.5-, 1-, 2-, and 3-liter solution bags for even and accurate heating. With the ivNow2, clinicians can easily customize setups by stacking or removing warming pods.

The unit installs on a countertop or mounts to the wall using the included mounting plate, providing flexibility for any workspace. Each chamber operates independently with its own electronic control and LED temperature display.

FEATURES

- Quickly warms and maintains temperature of IV, injection, and irrigation fluids
- Specially contoured warming module holds 0.5-, 1-, 2-, and 3-liter solution bags
- Stackable, field-customizable pod configuration for flexible setup
- Suitable for countertop or wall-mounted installation (mounting plate included)
- Individual electronic controls and LED display for each warming cavity

SPECIFICATIONS

- Model: P-2302-A
- Heating Time: Warms fluids to 104°F (40°C) in 30 minutes or less
- Fluid Capacity: Accommodates two (2) solution bags, each 0.5–3 liters
- Warming Output: Provides 52–120 liters of warmed fluid over a 12-hour period, depending on bag size
- Mounting Options: Suitable for countertop, wall-mounted, or pole-mounted installations

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)