

# Pedigo ivNow3 Fluid Warmer – New



Stock Photo – Options May Vary

## FACTORY NEW

Pedigo ivNow3 Fluid Warmer – New. The ivNow 3 Fluid Warmer, model P-2303-A, rapidly warms and maintains the temperature of intravenous, injection, and irrigation solutions before use. Its triple-chamber design lets clinicians warm three solution bags at once, improving workflow efficiency.

Each contoured warming cavity securely cradles 0.5-, 1-, 2-, and 3-liter solution bags for even and accurate heating every time. With flexible configuration options, ivNow3 can be stacked, expanded, or reduced to fit any clinical setting.

It easily installs on countertops or mounts to walls using the provided mounting plate. Each chamber features independent electronic controls and a bright LED temperature display, offering precise monitoring in either Celsius or Fahrenheit.

## FEATURES

- Rapidly warms and maintains the temperature of injection, intravenous, and irrigation solutions prior to use
- Specially contoured warming modules cradle 0.5-, 1-, 2-, and 3-liter solution bags
- Three-pod configuration warms up to three bags simultaneously
- Stackable, field-customizable design, pods can be added or removed as needed
- Flexible installation for countertop, wall, or boom mounting (mounting plate included)

## SPECIFICATIONS

- P-2303-A
- Heating Performance: Warms fluids to 104 °F (40 °C) in 30 minutes or less
- Capacity: Holds three (3) IV bags, 0.5 – 3 L each
- Warming Output: Delivers approximately 57 – 180 liters of warmed fluid in a 12-hour period
- Temperature Accuracy: Maintains within +0/-2 °C (+0/-3 °F) of set point

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)