

Blickman SS George Anesthesia Utility Table w/2 Drawers | 7754SS



Stock Photo - Options May Vary

FACTORY NEW

Blickman 7754SS - George Anesthesia Utility Table w/2 Drawers. A popular choice in the OR and throughout healthcare facilities as a general utility cart.

Blickman's superior construction adds a supporting channel underneath the top shelf to enhance its strength; a positive stop to prevent accidental drawer removal and ball-bearing slides for an easy-sliding drawer. The Blickman Ferguson Anesthesia Utility Table Travels easily on 3" casters.

FEATURES

- Popular in the OR and throughout healthcare facilities
- Used as a general utility cart
- Superior construction
- Stainless Steel
- Two Drawer

SPECIFICATIONS

- Model Number: 7754SS
- Part Number: 0227754000
- Two Drawers
- Stainless Steel Construction
- Dimensions: 20"L x 16"W x 34"H

**Brochure and images are a general overview of the product. Actual product options may vary.*