

# Pedigo Fluid & Blanket Warmer | P-2145



Stock Photo - Options May Vary

## FACTORY NEW

Pedigo Fluid & Blanket Warmer: P-2145. Pedigo leads in designing and constructing high-quality stainless steel and chrome products. Their products appear in emergency rooms, exam rooms, operating rooms, and sterile processing departments.

After more than a quarter-century, the Pedigo name remains synonymous with quality, durability, sterility, and excellent customer service. Pedigo Fluid & Blanket Warmer Features: 30 Liter total capacity. This combination of blanket and fluid warmers features a 4-digit LED display, a window door, and control lock programming.

It also includes auto time programming, a top compartment for fluids, and a lower compartment for blankets. The fluid compartment offers one easy-access epoxy-coated basket with a 30-liter bottle capacity.

## FEATURES

- Fluid and Blanket Warmer Combo
- 30-liter capacity
- One epoxy-coated basket with a maximum capacity of 30 liters
- Two selectable operating modes: 'INJECTION' or 'IRRIGATION'
- Separate temperature controls for each compartment

## SPECIFICATIONS

- Overall Size: 31" W x 27" D x 71" H
- Fluid Compartment Size: 30 liters
- Fluid Compartment Dimensions: 26 1/2"w x 23 1/4"d x 14 7/8"h
- Blanket Compartment Size: 11.875 cu. ft.
- Blanket Compartment Dimensions: 26 1/2"w x 23 1/4"d x 32 1/2"h

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free      Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)