

Blickman Saturn Anesthesia Utility Table w/3 Drawers | 7857SS



Stock Photo - Options May Vary

FACTORY NEW

Blickman 7857SS - SS Saturn Anesthesia Utility Table w/3 Drawers. A popular choice in the OR and throughout healthcare facilities as a general utility cart. Blickman's superior construction adds a supporting channel underneath the top shelf to enhance its strength; a positive stop to prevent accidental drawer removal and ball-bearing slides for an easy-sliding drawer.

The Blickman Saturn Anesthesia Utility Table travels easily on 3" casters. Choose one to four drawers depending on your requirements. Comes with a Lifetime Limited Warranty.

See below for more information.

FEATURES

- Anesthesia Utility Table
- Model number 7857SS
- Three Drawers
- Popular in the OR and throughout healthcare facilities
- General utility cart

SPECIFICATIONS

- Model Number: 7857SS
- Part Number: 0237857000
- Three Drawers
- Dimensions: 20"L x 16"W x 34"H
- Casters: 3" Wheel Diameter

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)