

Pedigo Fluid & Blanket Warmer | P-2148



Stock Photo - Options May Vary

FACTORY NEW

Pedigo P-2148 Fluid & Blanket Warmer. Pedigo leads in designing and constructing high-quality stainless steel and chrome products.

Their products are found in emergency rooms, exam rooms, operating rooms, and sterile processing departments. After more than a quarter of a century, the Pedigo name remains synonymous with quality, durability, sterility, and outstanding customer service.

Pedigo Fluid & Blanket Warmer: P-2148 Features: This combination of blanket & fluid warmers features a 4-digit LED display, a window door, control lock programming, and auto time programming—a top compartment for fluids, a lower compartment for blankets. The fluid compartment includes one easy-access, epoxy-coated basket with a maximum capacity of 18 1-liter bottles or equivalent.

FEATURES

- Fluid and Blanket Warmer Combo
- 18" (457mm) depth perfect for built-in OR applications
- One epoxy-coated basket with a maximum capacity of 18 liters
- Two selectable operating modes: 'INJECTION' or 'IRRIGATION'
- Separate temperature controls for each compartment

SPECIFICATIONS

- Overall Size: 31"W x 21¾"D x 66½"H
- Fluid Compartment Size: 29 liters
- Fluid Compartment Dimensions: 26 1/2"w x 17 3/4"d x 12 1/2"h
- Blanket Compartment Size: 8.7 cu. ft.
- Blanket Compartment Dimensions: 26 1/2"w x 17 3/16"d x 32 1/2"h

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)