

Pedigo Deluxe Warming Cabinet | P-2055

FACTORY NEW

Pedigo Deluxe Warming Cabinet: P-2055. Pedigo is a leader in designing and constructing high-quality stainless steel and chrome products.

Consequently, Pedigo products are found in emergency rooms, exam rooms, operating rooms, and sterile processing departments. After more than a quarter-century, the Pedigo name remains synonymous with quality, durability, sterility, and exceptional customer service.

Pedigo Deluxe Warming Cabinet Features: Deluxe combination blanket & fluid warmers feature a 4-digit LED display, a window door, control lock programming, and auto time programming. Pedigo Warmers Feature Radiant Heat Technology, a method for distributing heat balance throughout the cabinet, to ensure accurate set-point temperatures.



Stock Photo - Options May Vary

FEATURES

- Deluxe combination blanket & fluid warmers
- Radiant Heat Technology
- four digits LED display
- A windowed door
- Control lock programming and auto-time programming

SPECIFICATIONS

- Compartment Size: 20.6 cubic ft total
 - Compartment Dimensions: 26½" W x 23-3/4" D x 56-3/4" H
 - Overall Dimensions: 32" W x 29" D x 71" H
 - two doors, one compartment
 - 5" casters
- *Brochure and images are a general overview of the product. Actual product options may vary.

Toll-Free | Local
 800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)