

Blickman Solution Warming Cabinet | 7927TG



Stock Photo - Options May Vary

FACTORY NEW

Blickman Solution Warming Cabinet: 7927TG. The Blickman 7927TG Solution Warming Cabinet offers reliable performance in tight spaces. It features a compact design that fits easily on countertops and tabletops, and the 20 1/8-inch depth supports convenient placement in crowded clinical areas.

The glass door provides clear visibility, and the self-closing magnetic gasket creates a smooth, secure seal. Additionally, the adjustable tension door improves user control during frequent access.

You can add an optional adjustable shelf, and this flexibility helps you organize various solution containers with ease. The 7927TG includes the same high-quality control features found in larger Blickman warmers, so you maintain consistent performance without increasing your footprint.

FEATURES

- Single Compartment
- Counter Top Warming Cabinet
- Comes with a Standard Glass Door
- PID Digital controllers feature Tactile Membrane (Stainless Steel Dome)
- Full Numeric Keypads

SPECIFICATIONS

- Over-All Dimension: 24"W 20-5/8"D 24-1/2"H
- Amps: 7.5
- Wattage: 750
- Foot Capacity: 3.36
- Limited Lifetime Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)