

# Blickman Solution Warming Cabinet | 7927TS



Stock Photo - Options May Vary

## FACTORY NEW

Blickman Solution Warming Cabinet: 7927TS. The Blickman 7927TS Solution Warming Cabinet delivers dependable performance in tight clinical spaces. Its compact design supports facilities that need efficient warming without sacrificing valuable work area. The single-compartment layout simplifies organization, and the cabinet includes self-closing magnetic gaskets that create a secure and smooth seal.

Additionally, the adjustable tension door provides steady control during frequent access. You may add an optional adjustable shelf, and this feature increases flexibility for different container sizes. The cabinet measures 20 1/8 inches deep, so it fits easily on most countertops or tabletop surfaces.

Furthermore, the compact footprint improves workflow in areas that require maximum space efficiency.

## FEATURES

- Single Compartment
- Counter Top Warming Cabinet
- Stainless Steel Door
- PID Digital controllers feature Tactile Membrane (Stainless Steel Dome)
- Full Numeric Keypads

## SPECIFICATIONS

- Over-All Dimension: 24"W 20-5/8"D 24-1/2"H
- Amps: 7.5
- Wattage: 750
- Foot Capacity: 3.36
- Limited Lifetime Warranty

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)