

Blickman Warming Cabinet | 7925TG



Stock Photo - Options May Vary

FACTORY NEW

Blickman Warming Cabinet: 7925TG. The Blickman 7925TG delivers compact and reliable heating performance for medical facilities that need a dedicated solution for warming. This 24.5-inch unit includes a single glass door, and it fits easily on existing cabinetry or a Blickman mobile stand.

It warms solutions with consistent accuracy, and it offers the same quality and control features found in larger Blickman models. Additionally, the 7925TG includes an adjustable shelf, and it supports flexible storage for various container sizes. The full numeric keypad sits conveniently above the door, and it allows quick temperature adjustments.

Staff navigate settings easily, and they manage daily tasks with confidence.

FEATURES

- Single Compartment
- Single Glass Door
- PID Digital controllers feature Tactile Membrane (Stainless Steel Dome)
- Full Numeric Keypads
- Over-Temp audible and Visual Warnings

SPECIFICATIONS

- Overall Dimensions: 30"W x 26-1/2"D x 24-1/2"H
- Amps: 7.5
- Wattage: 750
- Foot Capacity (overall): 5.21
- Limited Lifetime Warranty

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)