

# Blickman Warming Cabinet | 7925TS



Stock Photo - Options May Vary

## FACTORY NEW

Blickman Warming Cabinet: 7925TS. The Blickman 7925TS Warming Cabinet delivers dependable performance for solution warming in medical and laboratory environments. It stands 24.5 inches high and fits easily on existing cabinets or a Blickman mobile stand.

The single stainless steel door adds durability, and the compact footprint supports efficient workspace planning. Additionally, the adjustable shelf offers flexible storage for a wide range of solutions.

This cabinet provides accurate temperature control, and the full numeric keypad simplifies operation. You access the keypad above the door, and you adjust settings quickly during fast-paced workflows.

## FEATURES

- Single Compartment
- Single Stainless Steel Door
- PID Digital controllers feature Tactile Membrane (Stainless Steel Dome)
- Full Numeric Keypads
- Over-Temp audible and Visual Warnings

## SPECIFICATIONS

- Overall Dimensions: 30" W x 26-1/2" D x 24-1/2" H
- Amps: 7.5
- Wattage: 750
- Foot Capacity (overall): 5.21
- Limited Lifetime Warranty

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)