

Blickman Solution/Blanket Warming Cabinet | 7922TG



Stock Photo - Options May Vary

FACTORY NEW

Blickman Solution/Blanket Warming Cabinet: 7922TG. delivers compact, high-performance heating for medical environments that need reliable warming in tight spaces. This 36-inch cabinet fits neatly under countertops, and it supports placement on existing cabinetry or a Blickman mobile stand.

It heats solutions or blankets with consistent accuracy, and it adds flexibility to clinics, surgery centers, and patient care areas. Additionally, the 7922TG includes a built-in temperature data recorder with a real-time clock, and the easy-to-use keypad sits directly above the door for quick access.

The graphical display features 160 x 180 pixels in dual colors, and it gives staff clear temperature visibility.

FEATURES

- Single Compartment and Single Glass Door
- PID Digital controllers feature Tactile Membrane (Stainless Steel Dome)
- Full Numeric Keypads
- Over-Temp audible and Visual Warnings & Automatic Cut-Off
- Built-In Temperature Data Recorder with Real-Time Clock

SPECIFICATIONS

- Model Number: 7922TG
- Part Number: 14B7922243
- Width: 30"
- Depth: 26.625"
- Height: 35.5"

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)