

# Blickman Solution/Blanket Warming Cabinet | 7922TS



Stock Photo - Options May Vary

## FACTORY NEW

Blickman 7922TS Solution/Blanket Warming Cabinet delivers dependable heating performance in a compact design. This 36-inch unit features a single stainless steel door and fits easily under countertops or in tight spaces.

It warms solutions or blankets with consistent accuracy, and the Blickman 7922TS supports fast access in clinics, surgery centers, and patient care areas. Moreover, the 7922TS includes a built-in temperature data recorder, and it offers flexible placement on existing cabinets or a Blickman mobile stand.

The keypad sits above the door, and it allows quick temperature adjustments. The self-tune feature maintains the set temperature within  $\pm 1^{\circ}\text{F}$ , ensuring reliable daily operation.

## FEATURES

- Single Compartment
- Single Stainless Steel Door
- PID Digital controllers feature Tactile Membrane (Stainless Steel Dome)
- Full Numeric Keypads
- Over-Temp audible and Visual Warnings

## SPECIFICATIONS

- Dimensions: 30"W 26-1/2"D 36"H
- Amps: 7.5
- Wattage: 750
- Foot Capacity: 9.04
- 2 Shelves

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)