

Blickman Stainless Steel Jumbo Mayo Stand | 8871SS



Stock Photo - Options May Vary

FACTORY NEW

Blickman 8871SS - Stainless Steel Jumbo Mayo Stand. Number one in ease and convenience. Just lift the tray frame - it automatically locks. To lower, hold the tray frame and press the foot-operated control button.

Release to lock in place. A guard over the foot button prevents accidental release. This Blickman Mayo Stand comes with a Limited Lifetime Warranty. Blickman's State of the Art Double Post Design minimizes lateral movement and increases overall unit stability.

Select our Model 8871SS, double-post jumbo Benjamin Mayo Stand when a larger tray size is required. The 8871SS is height adjustable from 39-1/2" to 62" with a tray size of 20" x 25".

FEATURES

- Just Lift the Tray Frame - It Automatically Locks
- To Lower - Hold Tray Fram and Press Foot Operated Control Button
- Release to Lock in Place
- Guard Over the Foot Button Prevents Accidental Release
- Stainless Steel Tray Included

SPECIFICATIONS

- Benjamin Stand
- Height Adjustment: 39.5" to 62"
- Tray: Stainless Steel
- Tray Dimensions: 20" x 25"
- Estimated Shipping Weight: 54 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)