

Blickman Stainless Steel Mayo Stand | 1520SS



Stock Photo - Options May Vary

FACTORY NEW

Blickman 1520SS - Stainless Steel Mayo Stand. These Newark stainless steel Mayo Stands are a superior choice for value-conscious customers.

They utilize a Tru-Loc™ friction knob for manual locking at the desired height. The 1520SS uses lighter tray support and a stainless steel base (one size only).

Additional Features: Blickman's Stainless Steel Mayo Stand model 1520SS travels easily on 2", single-wheel casters with a 37 1/4" - 64" height adjustment and comes with a Limited Lifetime Warranty. Note: Click the following link to view the wide variety of Mayo Stand models available to best meet your needs from Venture Medical.

FEATURES

- Superior Choice For Value-Conscious Customers
- Tru-Loc™ Friction Knob
- Lighter Tray Support
- Stainless Steel Base
- Includes Tray

SPECIFICATIONS

- Model Number: 1520SS
- Part Number: 0671520000
- Stainless Steel Base
- Single-Post Mayo Stand
- &

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)