

Blickman Stainless Steel Mayo Stand | 8867SS



Stock Photo - Options May Vary

FACTORY NEW

The Blickman 8867SS is a stainless steel Mayo stand designed for ease of use in clinical settings. The tray automatically locks when lifted and lowers via a foot-operated control button, with a guard preventing accidental release. The double-post design minimizes lateral movement and enhances overall unit stability.

This Benjamin-style stand features an adjustable height range and includes the tray. The stainless steel construction ensures durability and easy cleaning. For applications requiring a larger tray surface, the Model 8871SS double-post jumbo version is available.

FEATURES

- Just Lift the Tray, it Automatically Locks
- To Lower, Hold the Tray Frame and Press the Foot Operated Control Button
- Release to Lock in Place
- Guard Over the Foot Button Prevents Accidental Release
- Tray Included
- Benjamin Stand
- Adjustable Height

SPECIFICATIONS

- Height Adjustment Range: 39.5" to 62"
- Tray Construction: Stainless Steel
- Tray Dimensions: 12-5/8" x 19-1/8" tray

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)