

Blickman Stainless Steel Mayo Stand | 8869SS



Stock Photo - Options May Vary

FACTORY NEW

The Blickman 8869SS is a stainless steel Mayo stand featuring a foot-operated design for convenient height adjustment. Lifting the tray frame automatically locks it in place, while pressing the foot-operated control button allows for lowering. A guard over the foot button prevents accidental release during use.

The stand incorporates a state-of-the-art double post design that minimizes lateral movement and increases overall unit stability. This construction ensures reliable performance in clinical and surgical environments.

The unit comes with a stainless steel tray and is backed by a Limited Lifetime Warranty.

FEATURES

- Automatically Locks
- Guard Over the Foot Button Prevents Accidental Release
- Tray Included
- Adjustable Height
- Benjamin Stand

SPECIFICATIONS

- Height Adjustment Range; 39.5" to 62"
- Construction: Stainless Steel Tray
- Tray Dimensions: 16" x 21"
- Shipping Weight: 46 lbs.

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)