

# Blickman Stanford Stainless Steel Mayo Stand | 8848SS



Stock Photo - Options May Vary

## FACTORY NEW

Blickman 8848SS - Stanford Stainless Steel Mayo Stand. The 8848SS Mayo Stand is built with a heavy-duty stainless steel base and four casters for easy maneuverability, making it ideal for use in busy environments like ERs and exam rooms.

Its durable construction ensures long-lasting reliability, while the four caster wheels allow for smooth movement across various surfaces, even in tight spaces. The stand has a Tru-Loc™ friction knob for height adjustments from 39-1/2" to 62", offering flexibility for tasks.

The tray size is 12-5/8" x 19-1/8", providing ample space for instruments and supplies. This makes it a versatile and essential tool for healthcare professionals.

## FEATURES

- Easy Maneuverability
- Two Additional Casters
- Tubular Frame
- Ideal for ER or Exam Rooms
- Adjustable Height

## SPECIFICATIONS

- Height Adjustment Range: 39-1/2" to 62"
- Tray Dimensions: 12-5/8" x 19-1/8"
- Tray Construction: Stainless Steel
- Part Number: 0648848000
- Estimated Shipping Weight: 35 lbs.

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)