

Blickman Ultra Space Saver Case Cart | CCC5-19



Stock Photo - Options May Vary

FACTORY NEW

Blickman CCC5-19 - Ultra Space Saver Case Cart. This Ultra Space-Saver Case Cart offers greater capacity using the same floor space as the Multi-Purpose Cart.

It's useful for routine or extensive procedures requiring increased volume/size of instruments and/ or supplies. Body - 18 gauge, Door(s) - 18 gauge, Top - 16 gauge, Chassis - 11 gauge, Bumper channel - 14 gauge. Improved features include Extra Sound-Deadening Channels, all over Heavier Duty, 2 Push Handles, and 2 Additional Recessed Push Handles.

Additional Features: Stainless Steel welded, Angle Frame Chassis absorbs impact transmitted through wheels and bumper frames, thus protecting contents. Protects the body against distortion while preventing looseness and rattling.

FEATURES

- Ultra Space-Saver Case Carts
- Uses the same floor space as the Multi-Purpose Cart
- Ideal for routine or extensive procedures requiring increased volume/size of instruments and/ or supplies
- Stainless Steel welded, Angle Frame Chassis
- Removable pull-out type stainless steel wire shelves are provided on closed carts

SPECIFICATIONS

- Model Number: CCC5-19
- Part Number: 2293335000
- Overall Dimensions: 42"W x 29"D x 55-1/2"H
- Casters: 6" diameter s/s, 3 Swivel, 1 with Swivel Lock
- Shelves: 2 Interior

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)