

# Enthermics 6-Blanket Warming Cabinet | DC250



Stock Photo - Options May Vary

## FACTORY NEW

The Enthermics DC250 6-Blanket Warming Cabinet is an ideal solution for healthcare settings. Its countertop design fits seamlessly into labor and delivery rooms, ICUs, and NICUs. With no moving parts, the warmer operates quietly, allowing healthcare providers to focus on patient care without disruption.

This design also ensures ease of use, enabling nurses and aides to position it where needed to maximize efficiency. The DC250 incorporates features that enhance both comfort and convenience. Inside, relaxing blue lights create a soothing atmosphere for patients.

Additionally, WarmSafe™ zone heating technology ensures rapid and efficient warming of blankets. This technology, paired with heated shelves, provides consistent and even warmth throughout the cabinet.

## FEATURES

- LED interior lighting casts a comforting blue glow with two (2) different intensity settings and off-mode
- In the event of a power failure, the cabinet will remember its programming and begin to operate as before when power is restored
- A warming shut-off system, separate from the electronic control, prevents overheating
- The access point and the plug have been added to the back panel to allow the appliance owner to add their data logger or temperature management hardware
- Single-chamber warming cabinet White epoxy-coated steel exterior casing and interior insert

## SPECIFICATIONS

- Chamber: Single-chamber warming cabinet
- Lighting: LED
- Dimensions: 21.9 in H x 18.5 in W x 24.8 in D
- Interior Size: 2.5 ft
- Capacity: 4-6 Blankets

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free | Local  
800-627-3215 | 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)