

# Bovie Medical Hand-Held Light | UV59

## FACTORY NEW

Bovie Medical Hand-Held Light: UV59. Hand-held, lightweight, battery-operated, portable light designed for convenience and ease of use.



*Stock Photo - Options May Vary*

It features optional ultraviolet lighting, which makes it ideal for dermatological examinations and detecting skin conditions. The compact design allows for easy portability, ensuring that medical professionals can take it anywhere.

Whether used in the office or in the field, the UV59 provides reliable illumination and versatility for various medical applications. Call for pricing.

## FEATURES

- The UV Woods Light is a hand-held, lightweight, battery-operated, portable light with optional ultraviolet lighting.
- Used in diagnostic visualization and semen detection.
- This device is intended to be used by trained healthcare professionals.

## SPECIFICATIONS

- Fluorescent Lamp (replaceable ( UV59B): 365 nm
- LED: (light-emitting diode) Lamp (non-replaceable)
- Power: 6V (replaceable 4 AA Batteries)
- Operating Principle: Used in diagnostic visualization and semen detection.
- Patient Population: No restriction

*\*Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free 800-627-3215 | Local 813-740-8640

[salesteam@venturemedical.com](mailto:salesteam@venturemedical.com)

[View this product online →](#)