

Burton LED Medical UV Light with Magnifier | UV604



Stock Photo - Options May Vary

FACTORY NEW

The Burton UV604 LED Medical UV Light with Magnifier is a professional-grade inspection tool that combines UV and white LED lighting in one compact unit. Users can easily switch between light modes for improved visibility and accuracy during clinical examinations. The ergonomic design ensures comfort during extended use.

The device features a high-quality 5-diopter lens with 2.25x magnification, allowing professionals to examine fine details with exceptional clarity. The powerful LED system delivers consistent light while maintaining energy efficiency and long battery life. A rechargeable lithium-ion battery enables cordless operation for maximum flexibility.

The magnifier supports a wide range of applications, including dermatology, forensics, and industrial inspection. It provides reliable, long-term performance with consistent output and durable construction suitable for professional medical and diagnostic environments.

FEATURES

- Combines UV and white LED lighting in one device for versatile use
- Delivers consistent 365 nm UVA output for accurate inspections
- Enhances detail visibility with 2.25x magnification and 5-diopter lens
- Provides high color accuracy with 95 CRI white LED illumination
- Supports cordless operation with rechargeable lithium-ion battery
- Ensures long-lasting performance with up to 36,000-hour LED lifespan
- Improves workflow with lightweight, ergonomic handheld design
- Adapts to dermatology, forensic, medical, and industrial applications

SPECIFICATIONS

- Light Source: UV LED (365 nm UVA) + White LED
- Magnification: 2.25x (5-diopter lens)
- Illuminance: Up to 2,750–3,288 lux at 6 inches
- Color Temperature: 3800K (White LED)
- CRI: 95 (White LED)
- Viewing Area: 1.85" x 3.85"
- LED Lifespan: 36,000+ hours (white), 30,000+ hours (UV)
- Battery: Rechargeable Li-ion, ~4 hours runtime

**Brochure and images are a general overview of the product. Actual product options may vary.*

Toll-Free Local
800-627-3215 | 813-740-8640

salesteam@venturemedical.com

[View this product online →](#)